



comment	String
date	Date
dicomTime	Float
extension	String
fileType	String
frameNo	Int 32
height	Int 32
inDatabaseFolder	Bool
instanceNumber	Int 64
iPod	Bool
isKeyImage	Bool
modality	String
mountedVolume	Bool
name	String
numberOfFrames	Int 32
numberOfSeries	Int 32
path	String
rois	roi
series	Series
sliceLocation	Float
sopInstanceUID	String
stateText	String
type	String
width	Int 32

comment	String
date	Date
dateAdded	Date
dateOpened	Date
dicomTime	Float
displayStyle	Int 32
id	Int 32
images	Image
modality	String
mountedVolume	Bool
name	String
numberOfImages	Int 32
rotationAngle	Float
scale	Float
seriesDescription	String
seriesDICOMUID	String
seriesInstanceUID	String
seriesSOPClassUID	String
stateText	Int 32
study	Study
thumbnail	Binary
type	String
windowLevel	Float
windowWidth	Float
xFlipped	Bool
xOffset	Float
yFlipped	Bool
yOffset	Float

accessionNumber	String
albums	Album
comment	String
date	Date
dateAdded	Date
dateOfBirth	Date
dateOpened	Date
dicomTime	Float
dictateURL	String
expanded	Bool
hasDICOM	Bool
id	String
institutionName	String
modality	String
name	String
numberOfImages	Int 32
patientID	String
patientSex	String
patientUID	String
performingPhysician	String
referringPhysician	String
reportURL	String
series	Series
stateText	Int 32
studyInstanceUID	String
studyName	String
type	String