



Precision Diagnostic Systems, Inc.

QUOTATION NO: 2799 Version 02

DATE: 11/14/2010

This quotation is valid for 30 days

Precision Diagnostic Systems, Inc.

2N353 Landrose

Maple Park, IL 60151

Phone: (630) 365-1123 Fax: (630) 365-1124

[www.precisiondiagnostic.com](http://www.precisiondiagnostic.com)

E-mail: [rfnudo@precisiondiagnostic.com](mailto:rfnudo@precisiondiagnostic.com)

To:

Mitchell Coleman

Heritage Medical Imaging, Inc.

10836 Marietta Ave.

Culver City, CA 90232

[ninja2@ca.rr.com](mailto:ninja2@ca.rr.com)

Project:

Cancer Center of Irvine

Notes:

QUANT.	DESCRIPTION	UNIT PRICE	EXTENDED PRICE	NET AMOUNT
	<b>iQ-WEBX20 Server, DICOMREADER, iQ-VIEW Diagnostic Workstation, iQ-VIEW Administrative Workstation</b>			
1	<b><u>iQ-WEBX 20 Server with Integrated Web Distribution to include:</u></b>  DICOM 3.0 compliant PACS with fast DICOM archive and web-based viewing and managing interface; connection of to up to 20 DICOM 3.0 compliant imaging modalities and/or DICOM devices; no limitation in number of web users or number of concurrent users; Definable user privileges for radiologists, referring physicians, hospital staff, administrators; DICOM query/retrieve, DICOM print client, DICOM worklist client; configurable automatic image routing; no licensed limitation in number of DICOM images; MySQL database, Apache Web Server; Client support of Windows 2000/XP/7 computers with Internet Explorer 6 through 8; iQ ActiveX plug-in; DICOM Queriable Server <b>Server Hardware:</b>  Intel Quad Core Processor, 4GB RAM; Rack Mount Configuration 1.7 Terabyte Available Image Storage; RAID 5 Configuration, 5x1000GB Hot-swap drives Windows 2008 Server Operating System Integrated 10/100/1000 NIC Rapid Rails and cable management kit included			
1	<b><u>iQ-VIEW Diagnostic Workstation to include:</u></b>  Full FDA-Certified Workstation environment with full 3D image rendering, window/level, cross reference display, zoom/pan, scroll, cine. Patient CD Creation and importation of CD's from other institutions, DICOM Print to laser imagers, Paper Print to existing Windows-compatible paper printers. <b>Workstation Hardware:</b> Intel Core 2 Duo Processor, 2GB RAM 250 GB Storage Windows 7 Professional Operating System Dual 20" 3-Megapixel, High-Intensity, Portrait-Style, Greyscale Displays Integrated 10/100/1000 NIC, Extended Keyboard, Mouse, CD/DVD Burning Drive Uninterruptible Power Supply			
1	<b><u>iQ-VIEW Administrative Workstation to include:</u></b>  General-purpose Workstation for technologist and clinician image review with window/level, cross reference display, zoom/pan, scroll, cine. Patient CD Creation and importation of CD's from other institutions, DICOM Print to laser imagers, Paper Print to existing Windows-compatible paper printers. <b>Workstation Hardware:</b> Intel Core 2 Duo Processor, 2GB RAM 250 GB Storage Windows 7 Professional Operating System 24" Landscape DVI Display Integrated 10/100/1000 NIC, Extended Keyboard, Mouse, CD/DVD Burning Drive Imprinted CD's for Patient CD production, qty 100; Uninterruptible Power Supply			
1	<b><u>DICOMREADER</u></b>  Integrated in iQ-VIEW Administrative Workstation. Provides automated uploading of disks patients bring in from other institutions. Patient ID's may be changed at time of importation			
	<b><u>System Price</u></b>	<b>\$80,800.00</b>	<b>\$80,800.00</b>	<b>\$64,640.00</b>
1	<b><u>Delivery, installation and training</u></b>  All site preparation costs necessary to allow installation and use of these devices is the responsibility of the end user.  <b>Support &amp; Warranties:</b> <b>One year from first use, remote software support &amp; upgrades</b> <b>36-Month, same day, on-site hardware support from Dell on Server and Workstation computers</b> <b>Workstation 3-Megapixel Displays are warranted for a period of 36-Months parts from first use</b> <b>Normal (non-overnight) freight charges are included with quoted price</b>	<b>\$2,200.00</b>	<b>\$2,200.00</b>	<b>\$2,200.00</b>
	<b>Freight:</b>			
	Freight costs are included			
	<b>Terms:</b>	Sub-Total:	<b>\$83,000.00</b>	<b>\$66,840.00</b>
	50% downpayment with order/40% at shipment/10% upon completion	Handling:	-	-
	<b>Estimated Shipment:</b>	Sales Tax:	-	-
	04-06 weeks A.R.O.	<b>TOTAL PRICE:</b>	<b>\$83,000.00</b>	<b>\$66,840.00</b>